
MOLECULAR DIAGNOSTICS FOR SEXUAL TRANSMITTED DISEASES

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



Step 2 Paso 2

PCR*

* OpeGen kits include all the necessary reagents even Taq polymerase

Step 3 Paso 3

Hybridization protocol
Protocolo de hibridación

*25 µl Denaturing Solution + 25 µl PCR A+B in STD Panel Strip /
*25 µl Solución Desnaturalizante + 25 µl PCRA+B en STD Panel Strip

PRODUCTS

Productos

STD PANEL STRIP

Introduction

At present, sexually transmitted diseases (STDs) make up the most common group of infectious diseases that must be reported in the majority of the world's countries. Their incidence is elevated, and more than 20 pathogenic microorganisms are known to be transmitted by sexual contact, with *Chlamydia trachomatis*, *Neisseria gonorrhoeae*, *Trichomonas vaginalis* and *Mycoplasma (genitalium and hominis)* being the most common.

The current permanent trend of increasing sexually transmitted diseases, especially among young people under 25 years, makes their early and accurate diagnosis of great interest.

The STD Panel Strip test is a test based on the reverse blot technique that allows the detection of 10 pathogens associated to sexually transmitted diseases: *Chlamydia trachomatis* (discriminating variants L1, L2 and L3, that cause *lymphogranuloma venereum*), *Neisseria gonorrhoeae*, *Mycoplasma genitalium*, *Trichomonas vaginalis*, *Ureaplasma parvum*, *Ureaplasma urealyticum*, *Mycoplasma hominis*, *Herpes simplex 1*, *Herpes simplex 2*, and *Treponema pallidum*.

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
RESULTS

Resultados

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- Biernat-Sudolska. M., "The need to verify of positive *Mycoplasma hominis* results obtained using the *Mycoplasma* IST2 tests". *Journal of Laboratory Diagnostics* (2013), vol 49, pags 5-8.
- Domeika M. et al. "Guidelines for the laboratory diagnostics of *Chlamydia trachomatis* infections in East European countries". *JEADV* (2009). 23: 1353-1363.
- Sokolovskiy E. et al. "Guidelines for the laboratory diagnosis for syphilis in East European countries". *JEADV* (2010). 24: 1125-1134.
- Shipitsyna E. "Guidelines for the laboratory diagnosis of trichomoniasis in East European Countries". *Acta Derm. Venereol.* (2010). 90: 461-467.
- Savicheva A. et al. "Guidelines for laboratory diagnosis of *Neisseria gonorrhoeae* in East-European countries". *Acta Medica Lituanica* (2007). 14: 65-74.
- Domeika M. et al. "Guidelines for the laboratory diagnosis of genital herpes in eastern European countries". www.eurosurveillance.org.

Introduction

Human *Papillomavirus* (HPV) is responsible for a variety of skin and mucosa disorders. To date, over 100 different HPV genotypes have been described, of which more than 40 infect the anogenital tract. Of these, approximately 1/3 is associated with cervical cancer and anal neoplasm. The anogenital HPVs are generally divided into two categories: those with a low oncogenic potential risk (Low-Risk group) and those with a medium-high oncogenic potential risk (High-Risk group). The high-risk HPVs are generally associated with high-grade precancerous lesions and invasive cancer, while low-risk HPVs are frequently found in asymptomatic or benign conditions such as genital warts.

Intended use

High PapillomaStrip test allows qualitative detection of the DNA from 19 human *Papillomavirus* subtypes of medium-high risk: 16, 18, 26, 31, 33, 35, 39, 45, 51, 52, 53, 56, 58, 59, 66, 68, 69, 73 and 82.

On the other hand, Low PapillomaStrip test allows the genotyping of DNA samples from 18 human *Papillomavirus* subtypes of low risk: 6, 11, 40, 42, 43, 44, 54, 61, 62, 67, 70, 71, 72, 74, 81, 83, 84 y 91.

Specific primers for E6-E7 regions for each 37 HPVs detected with the kit are used to perform the multiplex PCR. Since the expression of the E6-E7 genes is directly associated with carcinogenesis, detection of these targets may provide significant biological and pathological information. Some studies showed that in cervical cancer samples that were positive for the E6 region, the L1 gene could not be detected in approximately 30% of cases, due to a gene deletion process that occurs when integrating DNA of high-risk HPV in epithelial cells.

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

PRODUCTS BASED ON OLIGOIMMUNOCHROMATOGRAPHY

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Colloid + 100 µl Hybridization B.
Coloide + 100 µl T. Hibridación

50 µl PCR + 50 µl Denaturing Sol.
+ 50 µl Sol. Desnaturalizante

Resuspended colloid
+ 10 µl denatured PCR

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

Productos

CT OligoGen

CT OligoGen is a test designed for the detection of *Chlamydia trachomatis* by simultaneous amplification of two independent targets of its DNA: one in the cryptic plasmid and another in its own genome and subsequent identification by immunochromatography. This test can detect infections caused by wild-type CT, the Swedish variant or those that do not carry the cryptic plasmid.

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
NG OligoGen

NG OligoGen is a test designed for the detection of *Neisseria gonorrhoeae* by simultaneous amplification of two independent targets of its DNA and subsequent identification in immunochromatographic strips.

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
Productos

MG OligoGen

MG OligoGen test is a test designed for the qualitative detection of *Mycoplasma genitalium* (MG) through the amplification of a specific genomic region of adhesin protein (MgPa) gene and its subsequent identification by immunochromatography. This test allows the simple, economical and reliable detection of all infections caused by MG.

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
TV OligoGen

TV OligoGen is a test designed for the qualitative detection of *Trichomonas vaginalis* (TV) through the amplification of a genomic region of TV (X83109), and its subsequent identification by immunochromatography. The test allows the simple, economical and reliable detection of all infections caused by the protozoan.

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1. Weinstock, H., S. Berman, and W. Cates Jr. 2004. Sexually transmitted diseases among American youth: incidence and prevalence estimates, 2000. *Perspect. Sex. Reprod. Health* 36(1): 6-10.
2. Soper, D. 2004. Trichomoniasis: under control or undercontrolled?. *American Journal of Obstetrics and Gynecology* 190 (1): 281-90.
3. Wolner-Hanssen, et al. 1989. Clinical manifestation of vaginal trichomoniasis. *JAMA*; 264:571-6.
4. Centers for Disease Control and Prevention. 2010. Sexually Transmitted Diseases Treatment Guidelines. *MMWR* 59: RR-12.
5. Wendel KA, Erbeling EJ, Gaydos CA, Rompalo AM. 2002. *Trichomonas vaginalis* polymerase chain reaction compared with standard diagnostic and therapeutic protocols for detection and treatment of vaginal trichomoniasis. *Clin Infect Dis.* 35(5): 576-80.

PRODUCTS

Productos


HSV 1/2

HSV OligoGen is a test designed for the qualitative detection and differentiation of *Herpes simplex virus* 1 and 2 through the simultaneous amplification of a common genomic region of DNA polymerase gene, and its subsequent identification by immunochromatography. The test allows for the simple, economical and reliable detection of all infections caused by the both types of virus.

Herpes simplex is an infectious viral disease, which is characterized by inflammation and the appearance of cutaneous lesions consisting of clusters of small vesicles surrounded by a red ring. These infections can affect the area around the mouth (oral herpes, also known as cold sores) or the area around the genitals (genital herpes). Until recently it was thought that oral herpes was caused by the *herpes simplex* type 1 virus and genital herpes by the *Herpes simplex* type 2 virus. Now we know that both viruses can cause lesions or sores in both parts of the body.

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1. Centers for Disease Control and Prevention. Sexually transmitted disease surveillance, 2008. Division of STD/HIV Prevention, Centers for Disease Control and Prevention, Atlanta, GA.
2. Ramaswamy M. et al, "Detection and typing of herpes simplex DNA in genital swabs by real time polymerase chain reaction", *Journal of Virological Methods* 126 (2005), 203-206.
3. Espy M. et al, "Diagnosis of herpes simplex virus infections in the clinical laboratory by LightCycler PCR", *Journal of Clinical Microbiology*, Vol 38, n° 2 (2000), 795-799.
4. Burrows J. et al, "Detection and subtyping of Hepres simplex virus in clinical samples by LightCycler PCR, enzyme immunoassay and cell culture", *BMC Microbiology* (2002), 2, 12.
5. Madhavan H. et al., "Detection of Herpes simplex virus (HSV) genome using polymerase chain reaction (PCR) in clinical samples. Comparison of PCR with standar laboratory methods for the deteccion of HSV", *Journal of Clinical Virology*, 14 (1999) , 145-151.

PRODUCTS

Productos

CT + LGV OligoGen

CT/LGV OligoGen is a test designed for the qualitative detection of *Chlamydia trachomatis* (CT) and their serotypes L1, L2 and L3, causing agents of Lymphogranuloma venereum (LGV), through the simultaneous amplification of two genomic regions, one in the cryptic plasmid and another in the genome of CT (omp gene), and the *pmpH* gen for LGV, and its subsequent identification by immunochromatography. The test allows for the simple, economical and reliable detection of all infections caused by the bacteria, whether or not (Swedish variant) it contains the cryptic plasmid.

Infections by *C. trachomatis* affect humans exclusively. Despite being the pathogen of the greatest prevalence in sexually transmitted infections (STIs) in the Western world, it is one of the least-diagnosed aetiologic agents, due to the high percentage of cases of asymptomatic infection.


There are 18 known serotypes distributed in three biovars or pathotypes, each responsible for a different infectious process. One of these three biovars is the Lymphogranuloma venereum (LGV): systemic disease that involves the invasion of the lymph nodes produced by serotypes L1, L2 and L3.

The treatment of Lymphogranuloma venereum (LGV) is longer than for other infections caused by *C. trachomatis* and consequently it is important to make a differentiation of this pathotype in order to make a proper treatment.

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1. ECDC Guidance- Chlamydia Control in Europe, Junio 2009.
2. Centers for Disease Control and Prevention. Sexually transmitted disease surveillance, 2008. Division of STD/HIV Prevention, Centers for Disease Control and Prevention, Atlanta, GA.
3. Ostergaard L. et al, "Microbiological aspect of the diagnosis of *Chlamydia trachomatis*", Best Practise & Research Clinical Obstetrics and Gynaecology, Vol 16, N° 6 pp 789-799, 2002.
4. Sstothard D.R. Et al, "Phylogenetic analysis of the *Chlamydia trachomatis* a major outer membrane protein and examination of potential pathogenic determinants", Infection and Immunity, Vol. 66 N° 8, pp 3618-3625, 1998.
5. Seth-Smith H.M.B. Et al, " Co-evolution of genomes and plasmids within *Chlamydia trachomatis* and the emergence in Sweden of a new variant strain", BMC Genomics, 10, 239, 2009.
6. Levett P. N. et al, "Evaluation of Three Automated Nucleic Acid Amplification Systems for Detection of *Chlamydia trachomatis* and *Neisseria gonorrhoeae* in First-Void Urine Specimens", Journal of Clinical Microbiology Vol 46 n°6 (2008), 2019-2111.
7. Parra-Sanchez M. et al, "Evaluation of the cobas 4800 CT/NG Test for detecting *Chlamydia trachomatis* and *Neisseria gonorrhoeae* DNA in urogenital swabs and urine specimens", Diagnostic Microbiology and Infectious Disease 74 (2012) 338-342.

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
CT/NG/MG OligoGen kit is a test designed for the detection of Chlamydia trachomatis (CT) Neisseria gonorrhoeae (NG) and Mycoplasma genitalium (MG), through the amplification of different genomic regions, specific for each pathogen, and its subsequent identification by immunochromatography. The test allows the simple, economical and reliable detection of all infections caused by CT, NG and MG.

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1. ECDC Guidance- Chlamydia Control in Europe, Junio 2009.
2. Centers for Disease Control and Prevention. Sexually transmitted disease surveillance, 2008. Division of STD/HIV Prevention, Centers for Disease Control and Prevention, Atlanta, GA.
3. Ostergaard L. et al, "Microbiological aspect of the diagnosis of Chlamydia trachomatis", Best Practise & Research Clinical Obstetrics and Gynaecology, Vol 16, N° 6 pp 789-799, 2002.
4. Shipitsyna E. "Guidelines for the laboratory diagnosis of Mycoplasma genitalium infections in East European Countries". Acta Derm. Venereol. (2010). 90: 461-467.
5. Sstothard D.R. Et al, "Phylogenetic analysis of the Chlamydia trachomatis a major outer membrane protein and examination of potential pathogenic determinants", Infection and Immunity, Vol. 66 N° 8, pp 3618-3625, 1998.
6. Seth-Smith H.M.B. Et al, " Co-evolution of genomes and plasmids within Chlamydia trachomatis and the emergence in Sweden of a new variant strain", BMC Genomics, 10, 239, 2009.
7. Liang Ma. "Genetic variation in the complete MgPa Operon and its repetitive chromosomal elements in clinical strains of Mycoplasma genitalium". PloS One (2010) Vol 5(12):e15660.
8. Goire N. et al, "A duplex Neisseria gonorrhoeae real time polymerase chain reaction assay targeting the gonococcal por A pseudogen and multicopy opa genes", Diagnostic Microbiology and Infectious Disease 61(2008), 6-12.
9. Jimmy Twin. "Comparison of two Mycoplasma genitalium real-time PCR detection methodologies" Journal of Clinical Microbiology (2011), Vol 49, N° 3 , 1140-1142.
10. Edger A. "A comparative study of three different PCR assays for detection of Mycoplasma genitalium in urogenital specimens from men and women" Journal of Medical Microbiology (2008) 57, 304-309.
11. Levett P. N. et al, "Evaluation of Three Automated Nucleic Acid Amplification Systems for Detection of Chlamydia trachomatis and Neisseria gonorrhoeae in First-Void Urine Specimens", Journal of Clinical Microbiology Vol 46 n°6 (2008), 2019-2111.
12. Chloé Le Roy. "The first performance report for the Bio-Rad Dx CT/NG/MG assay for simultaneous detection of Chlamydia trachomatis, Neisseria gonorrhoeae and Mycoplasma genitalium in urogenital samples". Journal of Microbiological Methods (2012) 89, 193-197.
13. Parra-Sanchez M. "Comparison of the CT OligoGen kit with cobas 4800 assay for detection of Chlamydia trachomatis" Enfermedades infecciosas y Microbiología clínica (2015) Vol 33, N° 10
14. Parra-Sanchez M. "Performance of the NG OligoGen kit for the diagnosis of Neisseria gonorrhoeae: comparison with cobas 4800 assay" Diagnostic Microbiology and Infectious Disease (2016), Vol.84, Issue 1, 4-6.

OPERON S. A.
Camino del Plano, 19
50410 Cuarte de Huerva (Zaragoza) • Spain
Tel: +34 976 503 597 • Fax: +34 976 503 531

e-mail address: sales@operon.es

www.operon.es