


INFECTIOUS DISEASE TESTING

State-of-the-Art Diagnosis

ELISA



Human

Diagnostics Worldwide

INFECTIOUS DISEASE ELISAs

Superior Quality for Manual or Automated Testing

Perfect for blood bank screening

Product	Cat. No.	Sensitivity	Specificity	Sample
HIV Ag/Ab	96 tests 51215	100 %	99,6 %	70 µl serum or plasma
new HIV Ag/Ab	480 tests* 512155	100 %	99,6 %	70 µl serum or plasma
anti-HIV-1/2	96 tests 51205C	100 %	99,6 %	30 µl serum or plasma
anti-HCV	96 tests 51275	100 %	99,7 %	100 µl serum or plasma
new anti-HCV*	96 tests 51250	100 %	99,7 %	100 µl serum or plasma
	480 tests 512505	100 %	99,7 %	100 µl serum or plasma
new HBsAg ultra sens*	96 tests 51050	100 %	99,6 %	100 µl serum or plasma
	480 tests 510505	100 %	99,6 %	100 µl serum or plasma
SYPHILIS SCREEN	96 tests 51005	100 %	100 %	50 µl serum or plasma

Reliable infection serology for prenatal care and routine testing

ToRCH	Cat. No.	Other Parameters	Cat. No.
CMV (Cytomegalo Virus)	IgG	EBV (Epstein-Barr Virus)	IgG
	IgM		IgM
HSV 1 (Herpes Simplex Virus 1)	IgG	H. PYLORI	IgG
HSV 2 (Herpes Simplex Virus 2)	IgG		IgA
HSV (Herpes Simplex Virus)	IgM	MEASLES VIRUS	IgG
RUBELLA VIRUS	IgG		IgM
	IgM µ-capture	MUMPS VIRUS	IgG
Toxo (Toxoplasma Gondii)	IgG		IgM
	IgM µ-capture	VZV (Varicella-Zoster Virus)	IgG
Hepatitis **			IgM
	anti-HBs	Antibody test	51255
anti-HBc	Antibody test	51260	CHLAMYDIA TRACHOMATIS
			IgG
			IgA
Tropical Diseases			
CHAGAS REC	IgG		51150
DENGUE	IgG		51240
	IgM		51140

Automation for all laboratories

- > The choice of three different instruments allows flexible automation for your laboratory.
- > All HUMAN assays pre-programmed!

* Available end of 2014
** Non CE marked products

